



Anti-IFNA1 (aa 24-133) polyclonal antibody (DPAB-DC1668)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The protein encoded by this gene is produced by macrophages and has antiviral activity. This gene is intronless and the encoded protein is secreted.
Immunogen	IFNA1 (NP_076918, 24 a.a. ~ 133 a.a) partial recombinant protein with GST tag. The sequence is CDLPETHSLDNRRITLMLLAQMSRISPSSCLMDRHDFGFPQEEFDGNQFQKAPASVLHEL IQQIFNLFTTKDSSAAWDEDLLDKFCTELYQQLNDLEACVMQEERVGETP
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	IFNA1 interferon, alpha 1 [Homo sapiens (human)]
Official Symbol	IFNA1

Synonyms	IFNA1; interferon, alpha 1; IFL; IFN; IFNA@; IFNA13; IFN-ALPHA; IFN-alphaD; interferon alpha-1/13; IeIF D; IFN-alpha 1b; IFN-alpha-1/13; interferon-alpha1; interferon alpha-D; interferon alpha 1b;
Entrez Gene ID	3439
Protein Refseq	NP_076918
UniProt ID	L0N195
Chromosome Location	9p22
Pathway	Autoimmune thyroid disease; Cytokine Signaling in Immune system; Cytokine-cytokine receptor interaction; Cytosolic DNA-sensing pathway