



Anti-SAX polyclonal antibody (DPAB-DC1704)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to fruit fly sax. The sequence is KKTTHKLSADEKIRLDFDEVC
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fruit fly
Conjugate	Unconjugated
Applications	WB, ICC, IP,
Format	Liquid
Size	100 µg
Buffer	In buffer containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	sax saxophone [Drosophila melanogaster (fruit fly)]
Official Symbol	SAX
Synonyms	SAX; saxophone; Bkr43E; Brk43E; CG1891; Dmel\CG1891; Sa; Sax; SAX; STK-B; CG1891-PA; CG1891-PB; CG1891-PC; sax-PA; sax-PB; sax-PC;

Entrez Gene ID	35731
Protein Refseq	NP_001246193
UniProt ID	Q7JQ36
Chromosome Location	43E18-43E18; 2-57 cM
Pathway	TGF-beta signaling pathway;
Function	ATP binding; protein binding; protein kinase activity; protein serine/threonine kinase activity