



Anti-SNF4 (native) polyclonal antibody (DPAB-DC1725)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Native purified Candida SNF4.
Isotype	IgG
Source/Host	Goat
Species Reactivity	Bacteria
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from 20 mM potassium phosphate buffer, 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C on dry atmosphere. After reconstitution with 0.1 mL of deionized water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	SNF4 release from glucose repression, invertase expression: allelic to possible fragments CaP19.5768 and CaP19.5769 [Candida albicans SC5314]
Official Symbol	SNF4

Synonyms	SNF4; hypothetical protein;
Entrez Gene ID	3640338
UniProt ID	Q5A8I0
Chromosome Location	chromosome: 6