



Anti-DNAJB13 (aa 71-180) polyclonal antibody (DPAB-DC1762)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	DNAJA13 belongs to the evolutionarily conserved DNAJ/HSP40 family of proteins. For background information on the DNAJ family, see 608375.[supplied by OMIM, Mar 2008]
Immunogen	DNAJB13 (NP_705842, 71 a.a. ~ 180 a.a) partial recombinant protein with GST tag. The sequence is GLKGGIPLEFGSQTPWTTGYVFHGKPEKVFHEFFGGNNPFSEFFDAEGSEVDLNFGLQG RGVKKQDPQVERDLYLSLEDLFFGCTKKIKISRRVLNEDGYSSTIKDKIL
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	DNAJB13 DnaJ (Hsp40) homolog, subfamily B, member 13 [Homo sapiens (human)]
Official Symbol	DNAJB13

Synonyms	DNAJB13; DnaJ (Hsp40) homolog, subfamily B, member 13; TSARG5; TSARG6; RSPH16A; dnaJ homolog subfamily B member 13; 1700014P03Rik; DnaJ-like protein; radial spoke 16 homolog A; DnaJ (Hsp40) related, subfamily B, member 13; testis and spermatogenesis cell-related protein 6; testis spermatogenesis apoptosis-related protein 6; testis spermatocyte apoptosis-related gene 6 protein; testis spermatogenesis apoptosis-related gene 3 protein; testis spermatogenesis apoptosis-related gene 6 protein;
Entrez Gene ID	374407
Protein Refseq	NP_705842
UniProt ID	P59910
Chromosome Location	11q13.4
Function	molecular_function; unfolded protein binding;