



Anti-USP17L2 (aa 431-530) polyclonal antibody (DPAB-DC1779)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	DUB3 is a member of the ubiquitin processing protease (UBP) subfamily of deubiquitinating enzymes. See USP1 (MIM 603478) for background information.[supplied by OMIM, Mar 2008]
Immunogen	DUB3 (NP_958804, 431 a.a. ~ 530 a.a) partial recombinant protein with GST tag. The sequence is ESTLDHWKFLQEQNKTKEFNVKVEGTLPDVLVIHQSKYKCGMKNHHPEQQSLLNLS STTPTHQESMNTGTLASLRGRARRSKGKNKHSKRALLVCQ
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	USP17L2 ubiquitin specific peptidase 17-like family member 2 [Homo sapiens (human)]
Official Symbol	USP17L2

Synonyms	USP17L2; ubiquitin specific peptidase 17-like family member 2; DUB3; DUB-3; USP17; ubiquitin carboxyl-terminal hydrolase 17; deubiquitinating enzyme 3; deubiquitinating protein 3; ubiquitin specific peptidase 17-like 2; ubiquitin thioesterase 17-like protein 2; deubiquitinating enzyme 17-like protein 2; ubiquitin thiolesterase 17-like protein 2; ubiquitin carboxyl-terminal hydrolase 17-like protein 2; ubiquitin-specific-processing protease 17-like protein 2;
Entrez Gene ID	377630
Protein Refseq	NP_958804
UniProt ID	Q6R6M4
Chromosome Location	8p23.1
Function	protein binding; ubiquitin thiolesterase activity; ubiquitin-specific protease activity;