



Anti-PSF1 (N-terminal) polyclonal antibody (DPAB-DC1789)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to N-terminus of Drosophila melanogaster Psf1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fruit fly
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Psf1 [Drosophila melanogaster (fruit fly)]
Official Symbol	PSF1
Synonyms	PSF1; CG9187; CG9187-PA; Dmel\CG9187; psf1; CG9187-PA; Psf1; Psf1-PA;

Entrez Gene ID	38136
Protein Refseq	NP_612086
UniProt ID	Q9W0I7
Chromosome Location	61F4-61F4
Pathway	Cell Cycle; DNA Replication; GINS complex; S Phase
Function	contributes_to 3-5 DNA helicase activity;