



Anti-JADE3 polyclonal antibody (DPAB-DC1793)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Specific to recombinant protein GX1028. This antibody detects mPHF16 protein.
Immunogen	Recombinant GST fusion protein corresponding to 162 amino acids of mouse Phf16.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	50 µg
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Jade3 jade family PHD finger 3 [Mus musculus (house mouse)]
Official Symbol	JADE3

Synonyms	JADE3; jade family PHD finger 3; Phf16; AI851988; AW544806; mKIAA0215; 5730598B06; protein Jade-3; PHD finger protein 16; jade family PHD finger protein 3;
Entrez Gene ID	382207
Protein Refseq	NP_001276613
UniProt ID	Q6IE82
Chromosome Location	X A1.3; X
Function	metal ion binding; molecular_function; zinc ion binding;