



# Anti-TLR12 (aa 700-800) polyclonal antibody (DPAB-DC1798)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	A synthetic peptide corresponding to amino acids 700-800 of mouse Tlr12.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB (Tissue lysate), Flow Cyt,
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	In PBS (0.05% BSA, 0.05% sodium azide)
<b>Preservative</b>	0.05% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">Tlr12 toll-like receptor 12 [ Mus musculus (house mouse) ]</a>
<b>Official Symbol</b>	TLR12
<b>Synonyms</b>	TLR12; toll-like receptor 12; Tlr11; Gm1365;

Entrez Gene ID	<a href="#">384059</a>
Protein Refseq	<a href="#">NP_991392</a>
UniProt ID	<a href="#">A2RRZ9</a>
Chromosome Location	4 D2.2; 4