



Anti-TLR12 (N-terminal) polyclonal antibody (DPAB-DC1799)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with KLH) corresponding to N-terminus mouse Tlr12.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA, Dot,
Format	Liquid
Size	100 µg
Buffer	In buffer containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Tlr12 toll-like receptor 12 [Mus musculus (house mouse)]
Official Symbol	TLR12
Synonyms	TLR12; toll-like receptor 12; Tlr11; Gm1365;

Entrez Gene ID	384059
Protein Refseq	NP_991392
UniProt ID	A2RRZ9
Chromosome Location	4 D2.2; 4