



Anti-SMSR polyclonal antibody (DPAB-DC1809)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to Drosophila SMSr.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fruit fly
Conjugate	Unconjugated
Applications	WB (Recombinant protein),
Format	Liquid
Size	100 µg
Buffer	In PBS (0.08% sodium azide)
Preservative	0.08% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	SMSr [Drosophila melanogaster (fruit fly)]
Official Symbol	SMSR
Synonyms	SMSR; CG32380; CG8572; CG8576; Css3beta; DmellCG32380; dmSMSr; dSMSr; CG32380-PA; CG32380-PB; CPE synthase; SMSr; SMSr-PA; SMSr-PB;

Entrez Gene ID	38823
Protein Refseq	NP_648114
UniProt ID	Q9VS60
Chromosome Location	65F7-65F9
Pathway	sphingomyelin metabolism.
Function	molecular_function;