



Anti-RNF207 (aa 135-244) polyclonal antibody (DPAB-DC1811)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	RNF207 (ring finger protein 207) is a protein-coding gene. GO annotations related to this gene include zinc ion binding.
Immunogen	FLJ46380 (NP_997279, 135 a.a. ~ 244 a.a) partial recombinant protein with GST tag. The sequence is HVFRGRPGSGFSSTSLGHLGPKCEPHYTGGETEVQNKGLEPVSQRWQRLRPFDLGRAHWS PIQGGVVDLHRRGSPVCRPGPTLKGLCYPSGIEAATAQGRWGQHAVPSGL
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	RNF207 ring finger protein 207 [Homo sapiens (human)]
Official Symbol	RNF207

Synonyms	RNF207; ring finger protein 207; C1orf188; RING finger protein 207;
Entrez Gene ID	388591
Protein Refseq	NP_997279
UniProt ID	Q6ZRF8
Chromosome Location	1p36.31
Function	zinc ion binding;