



Anti-RESP18 polyclonal antibody (DPAB-DC1814)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	RESP18 (regulated endocrine-specific protein 18) is a protein-coding gene. Diseases associated with RESP18 include pancreatic cancer, and pancreatitis.
Immunogen	Recombinant protein corresponding to amino acids of recombinant RESP18. The sequence is PTKTTESPLAKVNRDQCFTSEVSVKALKQEVANPVKITYRCSYGGLDMMQAPGPSKEEII YKIMRLLWATS
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC-P,
Format	Liquid
Size	100 µl
Buffer	In PBS, pH 7.2 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	RESP18 regulated endocrine-specific protein 18 [Homo sapiens (human)]
Official Symbol	RESP18
Synonyms	RESP18; regulated endocrine-specific protein 18; regulated endocrine-specific protein 18 homolog;
Entrez Gene ID	389075
Protein Refseq	NP_001007090
UniProt ID	Q5W5W9
Chromosome Location	2q35