



# Anti-XMRV XMVV62\_GP1 (N-terminal) polyclonal antibody (DPAB-DC1844)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	A synthetic peptide corresponding to N-terminus of XMRV XMVV62_gp1.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunoaffinity purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC-P,
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	In PBS, pH 7.6 (1% BSA, 0.1% sodium azide)
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Store at 4°C. Do not freeze.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">XMVV62_gp1 putative envelope polyprotein [ Xenotropic MuLV-related virus VP62 ]</a>
<b>Official Symbol</b>	XMVV62_GP1

<b>Synonyms</b>	XMVV62_GP1; putative envelope polyprotein;
<b>Entrez Gene ID</b>	<a href="#">3964425</a>
<b>Chromosome Location</b>	NC_007815.1