



Anti-MCCD1 polyclonal antibody (DPAB-DC1864)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	MCCD1 (mitochondrial coiled-coil domain 1) is a protein-coding gene. Diseases associated with MCCD1 include chronic venous leg ulcers, and systemic lupus erythematosus.
Immunogen	Recombinant protein corresponding to amino acids of recombinant MCCD1. The sequence is SQNPKASMEEQTNSRGNGKMTSPPRGPGTHRTAELARAEELLEQQLELYQALLEGQEGAW EAQALVLKIQKLKEQ
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC-P, IF,
Format	Liquid
Size	100 µl
Buffer	In PBS, pH 7.2 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	MCCD1 mitochondrial coiled-coil domain 1 [Homo sapiens (human)]
Official Symbol	MCCD1
Synonyms	MCCD1; mitochondrial coiled-coil domain 1; mitochondrial coiled-coil domain protein 1;
Entrez Gene ID	401250
Protein Refseq	NP_001011700
UniProt ID	P59942
Chromosome Location	6p21.33