



Anti-PARK (N-terminal) polyclonal antibody (DPAB-DC1870)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with KLH) corresponding to N-terminus of Drosophila park.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fruit fly
Purification	Ammonium sulfate precipitation
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µl
Buffer	In PBS (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	park parkin [Drosophila melanogaster (fruit fly)]
Official Symbol	PARK

Synonyms	PARK; parkin; CG10523; Dmel\CG10523; Dpark; dpk; SD01679; CG10523-PB; CG10523-PC; D-parkin; dparkin; park-PB; park-PC;
Entrez Gene ID	40336
Protein Refseq	NP_730600
UniProt ID	Q7KTX7
Chromosome Location	78C2-78C2
Pathway	Protein processing in endoplasmic reticulum; Ubiquitin mediated proteolysis;
Function	ubiquitin conjugating enzyme binding; ubiquitin-protein transferase activity; zinc ion binding;