



Anti-MDS1 (aa 80-169) polyclonal antibody (DPAB-DC1911)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	MDS1 (NP_004982, 80 a.a. ~ 169 a.a) partial recombinant protein with GST tag. The sequence is FELRESNMPGAGLGIWTKRKIEVGEKFGPYVGEQRSNLKDPSYGWEVHLPRSRVSVHSW LYLGKRSSDVGIAFSQADVMPGLQCAFLS
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	MDS1 myelodysplasia syndrome 1 [Homo sapiens (human)]
Official Symbol	MDS1
Synonyms	MDS1; myelodysplasia syndrome 1; PRDM3; MDS1-EVI1; myelodysplasia syndrome-associated protein 1;

Entrez Gene ID [4197](#)

Chromosome Location 3q26
