



Anti-CCAP (N-terminal) polyclonal antibody (DPAB-DC1926)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to N-terminus of fruit fly CCAP.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fruit fly
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Ccap Cardioacceleratory peptide [Drosophila melanogaster (fruit fly)]
Official Symbol	CCAP
Synonyms	CCAP; Cardioacceleratory peptide; ccap; CCAP; CG4910; DmCCAP; Dmel\CG4910; CCAP prepropeptide; CG4910-PA; CG4910-PB; Ccap-PA; Ccap-PB; cardioacceleratory peptide;

crustacean cardioactive peptide;

Entrez Gene ID	42683
Protein Refseq	NP_651083
UniProt ID	Q29QG1
Chromosome Location	94C1-94C1
Function	neuropeptide hormone activity; neuropeptide receptor binding;