



Anti-OS01G0674500 polyclonal antibody (DPAB-DC1937)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	It recognizes HSP in rice tissue lysate.
Immunogen	A synthetic peptide corresponding to Oryza sativa Os01g0674500.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rice
Conjugate	Unconjugated
Applications	WB (Tissue lysate),
Format	Liquid
Size	100 µl
Buffer	In serum
Preservative	None
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Os01g0674500 Os01g0674500 [Oryza sativa Japonica Group]
Official Symbol	OS01G0674500

Synonyms	OS01G0674500; hypothetical protein;
Entrez Gene ID	4325911
Protein Refseq	NP_001043840
UniProt ID	Q0JKH4
Chromosome Location	chromosome: 1