



# Anti-OS11G0163600 polyclonal antibody (DPAB-DC1942)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	It recognizes HSP in rice tissue lysate.
<b>Immunogen</b>	A synthetic peptide corresponding to Oryza sativa Os11g0163600.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Rice
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB (Tissue lysate),
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	In serum
<b>Preservative</b>	None
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">Os11g0163600 Os11g0163600 [ Oryza sativa Japonica Group ]</a>
<b>Official Symbol</b>	OS11G0163600

<b>Synonyms</b>	OS11G0163600; hypothetical protein;
<b>Entrez Gene ID</b>	<a href="#">4349867</a>
<b>Protein Refseq</b>	<a href="#">NP_001065832</a>
<b>UniProt ID</b>	<a href="#">Q53PP5</a>
<b>Chromosome Location</b>	chromosome: 11