



Anti-DJ-1BETA (N-terminal) polyclonal antibody (DPAB-DC1946)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with KLH) corresponding to N-terminus of Drosophila dj-1beta.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fruit fly
Purification	Ammonium sulfate precipitation
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µl
Buffer	In PBS (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	dj-1beta [Drosophila melanogaster (fruit fly)]
Official Symbol	DJ-1BETA

Synonyms	DJ-1BETA; CG1349; DJ-1; DJ-1b; DJ-1B; DJ-1beta; Dmel\CG1349; CG1349-PA; DJ-1B; dj-1beta; dj-1beta-PA;
Entrez Gene ID	43652
Protein Refseq	NP_651825
UniProt ID	Q9VA37
Chromosome Location	100A5-100A5