



Anti-RFPL4B (aa 127-226) polyclonal antibody (DPAB-DC1950)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	RFPL4B (ret finger protein-like 4B) is a protein-coding gene. GO annotations related to this gene include zinc ion binding. An important paralog of this gene is RFPL1.
Immunogen	LOC442247 (NP_001013756, 127 a.a. ~ 226 a.a) partial recombinant protein with GST tag. The sequence is TKDPLRLACVLGTPCFSSGQHYWEVEVGEVKSWSLGVCKEPADRKSNDLFPGHGFWISMKA GAIHANTHLERIPASPRRLRRVGIFLDADLEEIQFFDVDNN
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Cell lysate), WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	RFPL4B ret finger protein-like 4B [Homo sapiens (human)]
Official Symbol	RFPL4B

Synonyms	RFPL4B; ret finger protein-like 4B; RNF211; RING finger protein 211;
Entrez Gene ID	442247
Protein Refseq	NP_001013756
UniProt ID	Q6ZWI9
Chromosome Location	6q21
Function	zinc ion binding;