



Anti-NCAM2 (aa 598-695) polyclonal antibody (DPAB-DC2000)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The protein encoded by this gene belongs to the immunoglobulin superfamily. It is a type I membrane protein and may function in selective fasciculation and zone-to-zone projection of the primary olfactory axons.
Immunogen	NCAM2 (NP_004531, 598 a.a. ~ 695 a.a) partial recombinant protein with GST tag. The sequence is SIHGQPSSGKSFKLSITKQDDGGAPILEYIVKYRSKDKEDQWLEKKVQGNKDHILEHLQ WTMGYEVQITAANRLGYSEPTVYEFSPMPKPNIIKDTL
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [NCAM2 neural cell adhesion molecule 2 \[Homo sapiens \(human\) \]](#)

Official Symbol	NCAM2
Synonyms	NCAM2; neural cell adhesion molecule 2; NCAM21; NCAM-2; N-CAM-2;
Entrez Gene ID	4685
Protein Refseq	NP_004531
UniProt ID	O15394
Chromosome Location	21q21.1
Pathway	Cell adhesion molecules (CAMs);