



Anti-RPL26L1 (aa 1-100) polyclonal antibody (DPAB-DC2125)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a protein that shares high sequence similarity with ribosomal protein L26. It is not currently known whether the encoded protein is a functional ribosomal protein or whether it has evolved a function that is independent of the ribosome. Transcript variants utilizing alternative polyA signals exist.
Immunogen	RPL26L1 (NP_057177, 1 a.a. ~ 100 a.a) partial recombinant protein with GST tag. The sequence is MKFNPFVTSDRSKNRKRHFNAPSHVRRKIMSSPLSKELRQKYNVRSMPIRKDDEVQVVRG HYKGQQIGKVVQVYRKKYVIYIERVQREKANGTTVHVGIIH
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [RPL26L1 ribosomal protein L26-like 1 \[Homo sapiens \(human\) \]](#)

Official Symbol	RPL26L1
Synonyms	RPL26L1; ribosomal protein L26-like 1; RPL26P1; 60S ribosomal protein L26-like 1; ribosomal protein L26 homolog; ribosomal protein L26 pseudogene 1;
Entrez Gene ID	51121
Protein Refseq	NP_057177
UniProt ID	Q9UNX3
Chromosome Location	5q35.1
Pathway	Cap-dependent Translation Initiation; Eukaryotic Translation Elongation; Eukaryotic Translation Termination; GTP hydrolysis and joining of the 60S ribosomal subunit
Function	structural constituent of ribosome;