



Anti-RNFT1 (aa 1-100) polyclonal antibody (DPAB-DC2128)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	RNFT1 (ring finger protein, transmembrane 1) is a protein-coding gene. Diseases associated with RNFT1 include sertoli cell-only syndrome. GO annotations related to this gene include zinc ion binding. An important paralog of this gene is RNFT2.
Immunogen	LOC51136 (NP_057209, 1 a.a. ~ 100 a.a) partial recombinant protein with GST tag. The sequence is MQANRSQ LHSPPGTGSS EDASTPQCVHTRLTGEGSCPHSGDVHIQINSIPKECAENASSR NIRSGVHSCAHGCVHSRLRGHSHSEARLTDDTA AESGDHG
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [RNFT1 ring finger protein, transmembrane 1 \[Homo sapiens \(human\) \]](#)

Official Symbol	RNFT1
Synonyms	RNFT1; ring finger protein, transmembrane 1; PTD016; RING finger and transmembrane domain-containing protein 1;
Entrez Gene ID	51136
Protein Refseq	NP_057209
UniProt ID	Q5M7Z0
Chromosome Location	17q23.1
Function	zinc ion binding;