



Anti-RAPGEFL1 (aa 1-99) polyclonal antibody (DPAB-DC2135)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description RAPGEFL1 (Rap guanine nucleotide exchange factor (GEF)-like 1) is a protein-coding gene. GO annotations related to this gene include guanyl-nucleotide exchange factor activity. An important paralog of this gene is RAPGEF2.

Immunogen RAPGEFL1 (NP_057423, 1 a.a. ~ 99 a.a) partial recombinant protein with GST tag. The sequence is
MEGPEGLGRKQACLAMLLHFLDTYQGLLQEEEGAGHIIKDLYLLIMKDESLYQGLREDTL
RLHQLVETVELKIPEENQPPSKQVKPLFRHFRRIDSCLQ

Source/Host Mouse

Species Reactivity Human

Conjugate Unconjugated

Applications WB (Recombinant protein), ELISA,

Size 50 µl

Buffer 50 % glycerol

Preservative None

Storage Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [RAPGEFL1 Rap guanine nucleotide exchange factor \(GEF\)-like 1 \[Homo sapiens \(human\) \]](#)

Official Symbol	RAPGEFL1
Synonyms	RAPGEFL1; Rap guanine nucleotide exchange factor (GEF)-like 1; Link-GEFII; rap guanine nucleotide exchange factor-like 1; link GEFII; Link guanine nucleotide exchange factor II; RAP guanine-nucleotide-exchange factor (GEF)-like 1;
Entrez Gene ID	51195
Protein Refseq	NP_057423
UniProt ID	Q9UHV5
Chromosome Location	17q21.1
Function	guanyl-nucleotide exchange factor activity;