



Anti-ZNF219 (aa 143-221) polyclonal antibody (DPAB-DC2139)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ZNF219 (zinc finger protein 219) is a protein-coding gene. Diseases associated with ZNF219 include intrahepatic cholangiocarcinoma, and cholangiocarcinoma. GO annotations related to this gene include histamine receptor activity and sequence-specific DNA binding transcription factor activity. An important paralog of this gene is ZNF217.
Immunogen	ZNF219 (NP_057507, 143 a.a. ~ 221 a.a) partial recombinant protein with GST tag. The sequence is MQATPATEGLARPQAPSSSAFRCPYCKGKFRTSAERERHLHLHRPWKCGLCSFGSSQEE ELLHHSLTAHGAPERPLAA
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [ZNF219 zinc finger protein 219 \[Homo sapiens \(human\) \]](#)

Official Symbol	ZNF219
Synonyms	ZNF219; zinc finger protein 219; ZFP219;
Entrez Gene ID	51222
Protein Refseq	NP_001095142
UniProt ID	Q9P2Y4
Chromosome Location	14q11
Function	DNA binding; histamine receptor activity; metal ion binding; protein binding