



Anti-GLTP (aa 110-209) polyclonal antibody (DPAB-DC2140)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The protein encoded by this gene is similar to bovine and porcine proteins which accelerate transfer of certain glycosphingolipids and glyceroglycolipids between membranes. It is thought to be a cytoplasmic protein.
Immunogen	GLTP (NP_057517, 110 a.a. ~ 209 a.a) partial recombinant protein with GST tag. The sequence is SICDGERDENHPNLIRVNATKAYEMALKKYHGWIVQKIFQAALYAAPYKSDFLKALSKGQ NVTEEECLEKIRLFLVNYTATIDVIYEMYTQMNAELNYKV
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [GLTP glycolipid transfer protein \[Homo sapiens \(human\) \]](#)

Official Symbol	GLTP
Synonyms	GLTP; glycolipid transfer protein; truncated glycolipid transfer protein;
Entrez Gene ID	51228
Protein Refseq	NP_057517
UniProt ID	A0A024RBI7
Chromosome Location	12q24.11
Function	glycolipid binding; glycolipid transporter activity; lipid binding;