



Anti-PRR16 (full length) polyclonal antibody (DPAB-DC2157)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	PRR16 (proline rich 16) is a protein-coding gene. Diseases associated with PRR16 include smallpox, and endotheliitis.
Immunogen	LOC51334 (AAH38838.1, 1 a.a. ~ 234 a.a) full-length recombinant protein with GST tag. The sequence is MTDSSKTDLTNSSSSGTTASSLEKIKVQANAPLIKPPAHTSAILTVLRKPNPPPPPRLT PVKCEDPKRVVPTANPVKTNGTLLRNGGLPGGPNKIPNGDICCIPNSNLDKAPVQLLMHR PEKDRCPQAGPRERVRFNKVQYHGYCPDCDTRYNIKNREVHLHSEPVHPPGKIPHQGPP LPPTPHLPPFPLENGGMGISHSNSFPPIRPATVPPTAPKPQKTILRKSTTTTV
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Cell lysate), WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [PRR16 proline rich 16 \[Homo sapiens \(human\) \]](#)

Official Symbol	PRR16
Synonyms	PRR16; proline rich 16; DSC54; proline-rich protein 16; mesenchymal stem cell protein DSC54;
Entrez Gene ID	51334
Protein Refseq	NP_057728
UniProt ID	Q569H4
Chromosome Location	5q23.1