



Anti-ZFP644 polyclonal antibody (DPAB-DC2235)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Specific to recombinant protein GX0214. This antibody detects mZFP644 protein.
Immunogen	Recombinant GST fusion protein corresponding to 246 mouse Zfp644. The sequence is MMQNEEKYEKILKALNSRRIIPRFVAQKLSSGDDFLSHNVLP LDEYHNGLKTEALSVSA SEEEGLHFLSECGERKPELP SGRKNQSLT LIELLKSRRLGEERNSAVSPHKTHNQTARKR FVQKCVLPLNEDSPLIYQPQKMDLTMHSAIDCKQKKSRSRSGSKKKMLTLP HGADEVYIL RCRFCGLVFRGPLSVQEDWIKHLQRHIVNANLPRTGAGMVEVTSLLKKPASITETSF SLL MAEAAS
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB,
Format	Liquid
Size	50 µg
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Zfp644 zinc finger protein 644 [Mus musculus (house mouse)]
Official Symbol	ZFP644
Synonyms	ZFP644; zinc finger protein 644; Zep-2; BM-005; AA408882; D5Ertd689e; 1110068L01Rik; zinc finger motif enhancer binding protein 2;
Entrez Gene ID	52397
Protein Refseq	NP_081132
UniProt ID	E9Q624
Chromosome Location	5 E5; 5 51.43 cM
Function	molecular_function;