



Anti-PLAGL2 (aa 1-100) polyclonal antibody (DPAB-DC2265)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Pleiomorphic adenoma gene-like 2 is a zinc-finger protein that recognizes DNA and/or RNA.
Immunogen	PLAGL2 (NP_002648, 1 a.a. ~ 100 a.a) partial recombinant protein with GST tag. The sequence is MTTFFTSVPPWIQDAKQEEEVGWKLVPRPRGREAESQVKCQCEISGTPFSNGEKLPHSL PQPEQRPYSCPQLHCGKAFASKYKLYRHMATHSAQKPHQC
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Cell lysate), WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	PLAGL2 pleiomorphic adenoma gene-like 2 [Homo sapiens (human)]
Official Symbol	PLAGL2
Synonyms	PLAGL2; pleiomorphic adenoma gene-like 2; ZNF900; zinc finger protein PLAGL2; C2H2-type

zinc finger protein; pleiomorphic adenoma-like protein 2;

Entrez Gene ID	5326
Protein Refseq	NP_002648
UniProt ID	Q9UPG8
Chromosome Location	20q11.21
Function	RNA polymerase II core promoter proximal region sequence-specific DNA binding transcription factor activity involved in positive regulation of transcription; metal ion binding; sequence-specific DNA binding; sequence-specific DNA binding transcription fac