



Anti-PLGLB2 (aa 21-96) polyclonal antibody (DPAB-DC2275)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	PLGLB2 (plasminogen-like B2) is a protein-coding gene. Diseases associated with PLGLB2 include thrombosis. An important paralog of this gene is PLGLB1.
Immunogen	PLGLB2 (NP_002656, 21 a.a. ~ 96 a.a) partial recombinant protein with GST tag. The sequence is PLDDYVNTQGPSLFSVTKKQLGAGSREECAAKCEEDKEFTCRAFQYHSKEQQCVIMAENR KSSIIIRMRDAVLFEK
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	PLGLB2 plasminogen-like B2 [Homo sapiens (human)]
Official Symbol	PLGLB2

Synonyms	PLGLB2; plasminogen-like B2; PLGP1; plasminogen-like protein B; plasminogen pseudogene 1; type B plasminogen related; plasminogen-related protein B;
Entrez Gene ID	5342
Protein Refseq	NP_002656
UniProt ID	Q02325
Chromosome Location	2p11.2