



Anti-GPR84 (aa 208-316) polyclonal antibody (DPAB-DC2288)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	GPR84 (G protein-coupled receptor 84) is a protein-coding gene. Diseases associated with GPR84 include pertussis, and monocytic leukemia. GO annotations related to this gene include G-protein coupled receptor activity.
Immunogen	GPR84 (NP_065103, 208 a.a. ~ 316 a.a) partial recombinant protein with GST tag. The sequence is AAQALDQYKLRQASIHSNHVARTDEAMPGRFQELDSRLASGGPSEGISEPVSAATTQTL EGDSSEVGDQINSKRAKQMAEKSPPEASAKAQPIKGARRAPDSSEFGK
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [GPR84 G protein-coupled receptor 84 \[Homo sapiens \(human\) \]](#)

Official Symbol	GPR84
Synonyms	GPR84; G protein-coupled receptor 84; EX33; GPCR4; G-protein coupled receptor 84; inflammation-related G protein-coupled receptor EX33; inflammation-related G-protein coupled receptor EX33;
Entrez Gene ID	53831
Protein Refseq	NP_065103
UniProt ID	Q9NQS5
Chromosome Location	12q13.13
Pathway	GPCRs, Other.
Function	G-protein coupled receptor activity;