



Anti-CYTL1 (full length) polyclonal antibody (DPAB-DC2318)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C17 is a cytokine-like protein specifically expressed in bone marrow and cord blood mononuclear cells that bear the CD34 (MIM 142230) surface marker (Liu et al., 2000 [PubMed 10857752]).[supplied by OMIM, Mar 2008]
Immunogen	CYTL1 (AAH31391, 1 a.a. ~ 136 a.a) full-length recombinant protein with GST tag. The sequence is MRTPGPLPVLLLLLAGAPAAARPTPPTCYSRMRALSQEITRDFNLLQVSEPSEPCVRYLPR LYLDIHNYCVLDKLRDFVASPPCWKVAQVDSLKDKARKLYTIMNSFCRRDLVFLDDCNA LEYPIPVTTALPDRQR
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [CYTL1 cytokine-like 1 \[Homo sapiens \(human\) \]](#)

Official Symbol	CYTL1
Synonyms	CYTL1; cytokine-like 1; C17; C4orf4; cytokine-like protein 1; cytokine-like protein C17;
Entrez Gene ID	54360
Protein Refseq	NP_061129
UniProt ID	Q9NRR1
Chromosome Location	4p16-p15
Function	receptor binding;