



Anti-POU2AF1 (aa 159-255) polyclonal antibody (DPAB-DC2331)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	POU2AF1 (POU class 2 associating factor 1) is a protein-coding gene. Diseases associated with POU2AF1 include primary effusion lymphoma, and mediastinitis. GO annotations related to this gene include transcription coactivator activity and transcription cofactor activity.
Immunogen	POU2AF1 (NP_006226, 159 a.a. ~ 255 a.a) partial recombinant protein with GST tag. The sequence is PAVGPPLEGPEHQAPLTYFPWPQPLSTLPTSTLQYRPPAPALPGPQFVQLPISIPVQLQ DMEDPRRAASSLTIDKLLLEEDSDAYALNHTLSVEG
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [POU2AF1 POU class 2 associating factor 1 \[Homo sapiens \(human\) \]](#)

Official Symbol	POU2AF1
Synonyms	POU2AF1; POU class 2 associating factor 1; BOB1; OBF1; OCAB; OBF-1; POU domain class 2-associating factor 1; BOB-1; OCA-B; OCT-binding factor 1; B-cell-specific coactivator OBF-1; POU domain, class 2, associating factor 1;
Entrez Gene ID	5450
Protein Refseq	NP_006226
UniProt ID	Q16633
Chromosome Location	11q23.1
Function	DNA binding; protein binding; transcription coactivator activity; transcription cofactor activity