



Anti-RNF186 (aa 75-150) polyclonal antibody (DPAB-DC2338)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	RNF186 (ring finger protein 186) is a protein-coding gene. Diseases associated with RNF186 include ulcerative colitis. GO annotations related to this gene include zinc ion binding.
Immunogen	RNF186 (NP_061935, 75 a.a. ~ 150 a.a) partial recombinant protein with GST tag. The sequence is DNTWSITCPLCRKVTAVPGGLICSLRDHEAVVGQLAQPCTEVSLCPQGLVDPADLAAGHP SLVGEDGQDEVSAHV
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	RNF186 ring finger protein 186 [Homo sapiens (human)]
Official Symbol	RNF186

Synonyms	RNF186; ring finger protein 186; RING finger protein 186;
Entrez Gene ID	54546
Protein Refseq	NP_061935
UniProt ID	Q9NXI6
Chromosome Location	1p36.13
Function	zinc ion binding;