



Anti-GM-ASNASE1 (N-terminal) polyclonal antibody (DPAB-DC2359)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to N-terminus of soybean L-asparaginase.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Soybean
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	GM-ASNASE1 L-asparaginase [Glycine max (soybean)]
Official Symbol	GM-ASNASE1
Synonyms	GM-ASNASE1; L-asparaginase; L-asparagine amido hydrolase;

Entrez Gene ID	547703
Protein Refseq	NP_001236606
UniProt ID	Q8H2A9
Chromosome Location	chromosome: 14