



Anti-NUP62CL polyclonal antibody (DPAB-DC2366)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a protein containing a domain found in nucleoporins which are glycoproteins found in nuclear pore complexes. Alternative splicing results in multiple transcript variants.
Immunogen	Recombinant protein corresponding to amino acids of human NUP62CL. The sequence is VATPVMTYGHLEGLINEWNLELEDQEKYFLLQATQVNAWDHTLIENGEMIRILHGEVNVK KLDQKRLEQELDFILSQQLEFLTYLEESTRDQSGHLHYLQDADEEHVEISTRSAEF
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC-P, IF,
Format	Liquid
Size	100 µl
Buffer	In PBS, pH 7.2 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	NUP62CL nucleoporin 62kDa C-terminal like [Homo sapiens (human)]
Official Symbol	NUP62CL
Synonyms	NUP62CL; nucleoporin 62kDa C-terminal like; nucleoporin-62 C-terminal-like protein;
Entrez Gene ID	54830
Protein Refseq	NP_060151
UniProt ID	Q9H1M0
Chromosome Location	Xq22.3
Function	structural constituent of nuclear pore;