



Anti-CNTLN (aa 971-1070) polyclonal antibody (DPAB-DC2368)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description CNTLN (centlein, centrosomal protein) is a protein-coding gene. Diseases associated with CNTLN include seckel syndrome, and senior-loken syndrome. GO annotations related to this gene include phosphorelay sensor kinase activity.

Immunogen C9orf39 (NP_060208, 971 a.a. ~ 1070 a.a) partial recombinant protein with GST tag.
The sequence is
QNDVHVRRQIRELKKMKKNRDACKTSTHKAQTLAASILNISRSDL EEILD TEDQVEIEK
TKIDAENDKEWMLYIQKLLEGQSLALSPRLKCNGAIMAHQ

Source/Host Mouse

Species Reactivity Human

Conjugate Unconjugated

Applications WB (Recombinant protein), ELISA,

Size 50 µl

Buffer 50 % glycerol

Preservative None

Storage Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [CNTLN centlein, centrosomal protein \[Homo sapiens \(human\) \]](#)

Official Symbol	CNTLN
Synonyms	CNTLN; centlein, centrosomal protein; C9orf39; C9orf101; bA340N12.1; centlein; centrosomal protein;
Entrez Gene ID	54875
Protein Refseq	NP_001107867
UniProt ID	Q9NXG0
Chromosome Location	9p22.2
Function	molecular_function; phosphorelay sensor kinase activity;