



Anti-USP47 (aa 203-301) polyclonal antibody (DPAB-DC2381)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	USP47 (ubiquitin specific peptidase 47) is a protein-coding gene. GO annotations related to this gene include WD40-repeat domain binding and ubiquitin thiolesterase activity. An important paralog of this gene is USP24.
Immunogen	USP47 (NP_060414, 203 a.a. ~ 301 a.a) partial recombinant protein with GST tag. The sequence is TEQADLINELYQGKLDYVRCLECGYEGWRIDTYLDIPLVIRPYGSSQAFASVEEALHAF IQPEILDGPNQYFCERCKKKCDARKGLRFLHFPYLLTLQ
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [USP47 ubiquitin specific peptidase 47 \[Homo sapiens \(human\) \]](#)

Official Symbol	USP47
Synonyms	USP47; ubiquitin specific peptidase 47; TRFP; ubiquitin carboxyl-terminal hydrolase 47; ubiquitin thioesterase 47; deubiquitinating enzyme 47; ubiquitin thiolesterase 47; ubiquitin-specific-processing protease 47; Trf (TATA binding protein-related factor)-proximal homolog;
Entrez Gene ID	55031
Protein Refseq	NP_001269588
UniProt ID	Q68D97
Chromosome Location	11p15.3
Function	WD40-repeat domain binding; protein binding; ubiquitin-specific protease activity;