



Anti-ARMC4 (aa 945-1044) polyclonal antibody (DPAB-DC2389)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ARMC4 (armadillo repeat containing 4) is a protein-coding gene. Diseases associated with ARMC4 include alzheimer's disease. An important paralog of this gene is ARMC3.
Immunogen	ARMC4 (NP_060546, 945 a.a. ~ 1044 a.a) partial recombinant protein with GST tag. The sequence is AISRCCMWGRNRVAFGEHKAVAPLVRYLKSNDTNVHRATAQALYQLSEDADNCITMHENG AVKLLLDVMVGSPDQDLQEAAAGCISNIRRLALATEKARYT
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	ARMC4 armadillo repeat containing 4 [Homo sapiens (human)]
Official Symbol	ARMC4

Synonyms	ARMC4; armadillo repeat containing 4; CILD23; armadillo repeat-containing protein 4;
Entrez Gene ID	55130
Protein Refseq	NP_001276949
UniProt ID	B7Z7Y0
Chromosome Location	10p12.1-p11.23