



Anti-RFWD3 (aa 651-749) polyclonal antibody (DPAB-DC2390)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	RFWD3 (ring finger and WD repeat domain 3) is a protein-coding gene. GO annotations related to this gene include ubiquitin-protein ligase activity and MDM2/MDM4 family protein binding.
Immunogen	RFWD3 (AAH59371, 651 a.a. ~ 749 a.a) partial recombinant protein with GST tag. The sequence is DTGNPICSCQPVHTFFGGPTCKLLTKNAIFQSPENDGNILVCTGDEAANSALLWDAASGS LLQLDQTDQPVLDCPFEVNRNSYLATLTEKMHVHTYKWE
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	RFWD3 ring finger and WD repeat domain 3 [Homo sapiens (human)]
Official Symbol	RFWD3

Synonyms	RFWD3; ring finger and WD repeat domain 3; RNF201; E3 ubiquitin-protein ligase RFWD3; RING finger protein 201; RING finger and WD repeat domain-containing protein 3;
Entrez Gene ID	55159
Protein Refseq	NP_060594
UniProt ID	Q6PCD5
Chromosome Location	16q23.1
Function	MDM2/MDM4 family protein binding; ligase activity; p53 binding; protein binding