



# Anti-ILK-2 (C-terminal) polyclonal antibody (DPAB-DC2410)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	A synthetic peptide (conjugated with KLH) corresponding to C-terminus of human ILK-2.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Purification</b>	Protein G purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB (Cell lysate), ELISA,
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	In PBS (0.09% sodium azide)
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">ILK-2 integrin-linked kinase-2 [ Homo sapiens (human) ]</a>
<b>Official Symbol</b>	ILK-2

**Synonyms** ILK-2; integrin-linked kinase-2;

---

**Entrez Gene ID** [55522](#)

---