



Anti-DEPDC1 (aa 93-186) polyclonal antibody (DPAB-DC2421)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	DEPDC1 (DEP domain containing 1) is a protein-coding gene. GO annotations related to this gene include protein binding and GTPase activator activity. An important paralog of this gene is DEPDC4.
Immunogen	DEPDC1 (NP_060249, 93 a.a. ~ 186 a.a) partial recombinant protein with GST tag. The sequence is GSENVDDNNQLFRFPATSPKTLPRRYPELRKNNIENFSKDKDSIFKLRNLSRRTPKRHG LHLSQENGEEKIKHEIINEDQENAIIDNRELSQEDV
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [DEPDC1 DEP domain containing 1 \[Homo sapiens \(human\) \]](#)

Official Symbol	DEPDC1
Synonyms	DEPDC1; DEP domain containing 1; DEP.8; SDP35; DEPDC1A; DEPDC1-V2; DEP domain-containing protein 1A; cell cycle control protein SDP35;
Entrez Gene ID	55635
Protein Refseq	NP_001107592
UniProt ID	Q5TB30
Chromosome Location	1p31.2
Function	GTPase activator activity; protein binding;