



Anti-TEX11 (internal region) polyclonal antibody (DPAB-DC2502)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene is X-linked and is expressed in only male germ cells. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
Immunogen	A synthetic peptide corresponding to amino acids at internal region of human TEX11. The sequence is HNFLKQTGTFSNDSC
Source/Host	Goat
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	ELISA,
Format	Liquid
Concentration	0.5 mg/mL
Size	100 µg
Buffer	In 0.5 mg/mL Tris saline, pH 7.3 (0.02% sodium azide, 0.5% BSA)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	TEX11 testis expressed 11 [Homo sapiens (human)]
Official Symbol	TEX11
Synonyms	TEX11; testis expressed 11; TGC1; TSGA3; testis-expressed sequence 11 protein; testis expressed sequence 11;
Entrez Gene ID	56159
Protein Refseq	NP_001003811
UniProt ID	Q8IYF3
Chromosome Location	Xq13.1
Function	protein binding;