



Anti-FMN2 (aa 144-243) polyclonal antibody (DPAB-DC2531)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Formin homology (FH) domain proteins (see FMN1; MIM 136535) play a role in cytoskeletal organization and/or establishment of cell polarity.[supplied by OMIM, Apr 2010]
Immunogen	FMN2 (AAH14364.2, 144 a.a. ~ 243 a.a) partial recombinant protein with GST tag. The sequence is DQEAENSLETTHKCFLETTAYFFMKPKLGEKEVSPNAFFSIWHEFSSDFKDFWKKENKL LLQERVKEAEEVCRQKKGKSLYKIKPRHDSGIKAKISMKT
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Cell lysate), WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	FMN2 formin 2 [Homo sapiens (human)]
Official Symbol	FMN2

Synonyms	FMN2; formin 2; formin-2;
Entrez Gene ID	56776
Protein Refseq	NP_064450
UniProt ID	Q9HBL1
Chromosome Location	1q43
Pathway	Dorso-ventral axis formation;
Function	actin binding; molecular_function;