



Anti-NCLN (aa 44-144) polyclonal antibody (DPAB-DC2541)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	NCLN (nicalin) is a protein-coding gene. GO annotations related to this gene include protein binding and peptidase activity.
Immunogen	NCLN (NP_064555, 44 a.a. ~ 144 a.a) partial recombinant protein with GST tag. The sequence is HEFTVYRMQQYDLQGQPYGTRNAVLNTEARTMAAEVLSRRCVLMRLLDfsyeQYQKALRQ SAGAVVILPRMAAVPQDVVRQFMEIEPEMLAMETAVPVY
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Cell lysate), WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	NCLN nicalin [Homo sapiens (human)]
Official Symbol	NCLN

Synonyms	NCLN; nicalin; NET59; nicastrin-like protein; nicalin homolog (zebrafish);
Entrez Gene ID	56926
Protein Refseq	NP_064555
UniProt ID	Q969V3
Chromosome Location	19p13.3
Function	peptidase activity; protein binding;