



Anti-C11ORF30 (aa 1081-1178) polyclonal antibody (DPAB-DC2544)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C11orf30 (chromosome 11 open reading frame 30) is a protein-coding gene. Diseases associated with C11orf30 include sporadic breast cancer, and spondylitis.
Immunogen	C11orf30 (NP_064578, 1081 a.a. ~ 1178 a.a) partial recombinant protein with GST tag. The sequence is PASSEKQTASQVEQPIITQGSSVTKITFEGRPPTVTKITGGSSVPKLTSPVTSISPIQA SEKTAVSDILKMSLMEAQIDTNVEHMIVDPPKKALATS
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	C11orf30 chromosome 11 open reading frame 30 [Homo sapiens (human)]
Official Symbol	C11ORF30

Synonyms	C11ORF30; chromosome 11 open reading frame 30; EMSY; GL002; protein EMSY;
Entrez Gene ID	56946
Protein Refseq	NP_064578
UniProt ID	Q7Z589
Chromosome Location	11q13.5