



Anti-TWSG1 (aa 124-223) polyclonal antibody (DPAB-DC2560)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TWSG1 (twisted gastrulation BMP signaling modulator 1) is a protein-coding gene. Diseases associated with TWSG1 include holoprosencephaly.
Immunogen	TWSG1 (NP_065699, 124 a.a. ~ 223 a.a) partial recombinant protein with GST tag. The sequence is AEELSHHENLVSFLETVNQPHHQNVSVPSNNVHAPYSSDKEHMCTVVYFDDCMSIHQCKI SCESMGASKYRWFHNACCECIGPECIDYGSKTVKCMNCMF
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	TWSG1 twisted gastrulation BMP signaling modulator 1 [Homo sapiens (human)]
Official Symbol	TWSG1

Synonyms	TWSG1; twisted gastrulation BMP signaling modulator 1; TSG; twisted gastrulation protein homolog 1; twisted gastrulation homolog 1;
Entrez Gene ID	57045
Protein Refseq	NP_065699
UniProt ID	Q9GZX9
Chromosome Location	18p11.3