



Anti-REXO4 (aa 323-422) polyclonal antibody (DPAB-DC2566)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	REXO4 (REX4, RNA exonuclease 4 homolog (<i>S. cerevisiae</i>)) is a protein-coding gene. GO annotations related to this gene include exonuclease activity and sequence-specific DNA binding transcription factor activity. An important paralog of this gene is ISG20.
Immunogen	XPMC2H (NP_065118, 323 a.a. ~ 422 a.a) partial recombinant protein with GST tag. The sequence is LVGHALHNDLKVLFLDHPKKKIRDTQKYKPFKSQVKSGRPSLRLLSEKILGLQVQQAHC SIQDAQAAMRLYVMVKKEWESMARDRRPLLTAPDHCSDDA
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [REXO4 REX4. RNA exonuclease 4 homolog \(*S. cerevisiae*\) \[*Homo sapiens* \(human\) \]](#)

Official Symbol	REXO4
Synonyms	REXO4; REX4, RNA exonuclease 4 homolog (S. cerevisiae); REX4; XPMC2; XPMC2H; RNA exonuclease 4; exonuclease XPMC2; XPMC2 prevents mitotic catastrophe 2 homolog; Xenopus prevents mitotic catastrophe 2 homolog;
Entrez Gene ID	57109
Protein Refseq	NP_001266278
UniProt ID	B4E331
Chromosome Location	9q34.2
Function	exonuclease activity; poly(A) RNA binding; sequence-specific DNA binding transcription factor activity;