



Anti-SALL4 (aa 954-1053) polyclonal antibody (DPAB-DC2570)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The protein encoded by this gene may be a zinc finger transcription factor. Defects in this gene are a cause of Duane-radial ray syndrome (DRRS).
Immunogen	SALL4 (NP_065169, 954 a.a. ~ 1053 a.a) partial recombinant protein with GST tag. The sequence is PKEILAPSVNVDPVVWNQYTSMLNGLAVKTNEISVIQSGGVPTLPVSLGATSVVNNATV SKMDGSQSGISADVEKPSATDGVPKHQFPHFLEENKIAVS
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	SALL4 spalt-like transcription factor 4 [Homo sapiens (human)]
Official Symbol	SALL4

Synonyms	SALL4; spalt-like transcription factor 4; DRRS; HSAL4; ZNF797; dJ1112F19.1; sal-like protein 4; zinc finger protein 797; zinc finger protein SALL4;
Entrez Gene ID	57167
Protein Refseq	NP_065169
UniProt ID	Q9UJQ4
Chromosome Location	20q13.2
Pathway	Developmental Biology; Transcriptional regulation of pluripotent stem cells.
Function	DNA binding; metal ion binding;